




Desert Radiology intelConnect EV Online Image Viewer

Quick Reference Guide




Launch the DR intelConnectEV online image viewer via link:


<https://pacsportal.desertrad.com>

Login with username and password provided to you


If you do not have credentials, please contact Marketing for access at: 702-759-7775 or marketing@desertrad.com



- After login, the PATIENT dashboard (below) will appear, showing a list of provider patients referred in the last 7 days.
- You can see a larger date range of patients by selecting from the drop down list titled DATE RANGE.
- Please note: if a clinical team member (MA, scheduler, etc.) NOT PROVIDER is logged in this screen will appear empty.
- *Clinical team members will need to search for patients by name (see below image for more instruction).*



- To search by patient, enter the patient's name in the top left search bar (highlighted in yellow).
- You can search by the first two letters of the patient's last or first name and DOB. Search by DOB ONLY is not available. *Please note: If the patient has a compound surname such as De Garcia, type the last name followed by comma.*
- To view reports/images of patients referred to you by another provider for continuum of care, you **MUST** select "Include patients in which access is restricted." Please note, this box will not remain checked and needs to be selected for each patient search.



Quick Reference Guide

- To open a patient "chart" simply select and click on the exam study line that you wish to open (image A)
- Once selected, the patient chart (see image B) will open providing access to both the patient report and images for the selected exam.

Image A

The screenshot shows the 'My patients' section of the Desert Radiology intelConnect EV interface. At the top, there is a search bar and a welcome message for 'JOHN DOCTOR MD'. Below the search bar, a table lists 8 patients found. The table has columns for patient name, DOB, gender, exam type, and date. The first patient is 'PATIENT, ABC' with DOB 10/25/1983, and the second is 'PATIENT, DEF' with DOB 05/23/1972. The table also includes a 'Date range' dropdown set to 'Last 7 days' and a 'Sort by' dropdown set to 'Patient name (A-Z)'. A 'Refresh' button is located at the bottom right of the table.

Image B

The screenshot shows the 'Patient, ABC' chart in the Desert Radiology intelConnect EV interface. The patient's name, DOB (10/25/1983), gender (Female), and patient ID are displayed at the top. Below this, there is a 'View PDF...' button and a section for the patient's report. The report is titled 'CR XR C SPINE W FLEX EXT' and includes details such as 'Priority: Routine', 'Accession Number: 510194912', 'Referring Physician: NELLIS', 'Status: Final Report (ZZ)', 'Study Date: 03/28/2018 08:50 AM', 'Radiologist: GABAEFF MD, DINA', and 'Location: NELLIS'. The report text describes the examination as 'XR CERVICAL SPINE WITH FLEXION EXTENSION' and provides a 'TECHNIQUE' section. To the right of the report, there is a grid of 10 thumbnail images showing various views of the cervical spine, including AP, lateral, and oblique views.

- From the "patient chart" screen you can view and print a copy of the patient's report.
- To print the report, simply select VIEW PDF and follow the steps to print.
- To view the images with advanced tools, just click on the thumbnail on this screen to launch the image viewer.

